

Specifications CX-1

Dimensions	320 W x 531 D x 577 H mm, 26 kg	Focus Adjustment	Split Lines
Angle of view	MYD: 50 degrees, Non-Myd: 45 degrees 2 X magnification (digital)	Working distance	Corneal Reflection dots adjustment
Minimum pupil size	Myd: \varnothing 5.1mm (SP mode \varnothing 4.3 mm) Non-Myd: \varnothing 4.3 mm (SP mode \varnothing 3.8 mm)	Panning and Tilting	30 degrees to the left and right tilting range 15 degrees up, 10 degrees down
Working Distance	35mm	Light sources	Xenon tube for photography Halogenlamp for observation (Myd mode) IRED LED for observation (Non-Myd mode)
Photography modes	Colour /FA /Red Free /Cobalt and FAF	Fixation targets	External Internal LED dot matrix for Non-Myd mode (70 points) Internal fixation target for Myd mode (optional)
Mounted camera	Dedicated digital EOS camera (24 MegaPixel for current model) HDMI Output for external monitor Full HD resolution	Optional accessories	Stereo Unit SU-1 Internal eye fixation (CX-IF) Chin rest paper (500 sheets)
Flash levels	168 steps		
Patient's diopter	-31D ~ -7D, -10D ~+15D (standard) compensation +11D ~+33D		

Your Canon dealer:

Canon

CANON MEDICAL SYSTEMS EUROPE B.V.

<https://eu.medical.canon>

©Canon Medical Systems Corporation 2019. All rights reserved.
Design and specifications subject to change without notice.
Model number: CX1
MCAEC0002EUC 2019-08 CMSE/Printed in Europe

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485.
Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Made for Life is a trademark of Canon Medical Systems Corporation.

Disclaimer: Some features presented in this brochure may not be commercially available on all systems shown or may require the purchase of additional options. Please contact your local Canon Medical Systems representative for details.

Made For life